
# Puregraft Instrument Set 550/IS (SST)

- Prepare a fresh batch of the detergent according to the detergent manufacturer's recommendations.
- Fully immerse the devices into an ultrasonic bath and sonicate for a minimum of 10 min.-utes.
- Thoroughly rinse the devices with Reverse Osmosis (RO) or Deionized (DI) water to remove all detergent residues. Use a syringe to aid in rinsing.
- Dry the device(s) with a lint free cloth. Filtered, pressurized air can be used if applicable.
- Visually inspect each device for any visible soil.

## 2. STERILIZATION

The Puregraft Instrument Set (SST) components are provided non-sterile. Sterilization of the Puregraft Instrument Set (SST) components is required prior to use. For steam autoclaves, the Manufacturer has verified the following pre-vacuum sterilization cycles:

STERILIZER CYCLE TYPE	PRE-VACUUM
Temperature	min. 132 °C
Cycle Time	min. 3 minutes
Minimum Dry Time	20 minutes
Sample Configuration	Double Wrapped

It is the facility's responsibility to verify the appropriateness of the sterilization cycle, taking into consideration bioburden, microorganism resistance to heat (D-value), equipment loading configurations and heat distribution, the appropriate sterility assurance level, and other sterilization cycle parameters. The Manufacturer recommends that the facility establishes its own sterilization cycles by validating the sterilization equipment and cycles through use of appropriate standards. Do not exceed 138°C cycle temperature or 20 minute cycle time.

## INSTALLATION

Assemble easel by placing easel stand on flat surface. Install Easel plate onto Easel stand as shown until Easel plate is seated properly on stand as shown on IFU cover page.

## STORAGE AND HANDLING

- Protect from moisture during storage.

## WARNINGS AND PRECAUTIONS

- The physician should be familiar with operation of the device(s) and surgical procedure prior to performing the surgery.
- The Puregraft Instrument Set (SST) is provided non-sterile and MUST be cleaned and sterilized prior to use.
  - Inspect for cleanliness and damage prior to use.
  - DO NOT use the Puregraft Instrument Set (SST) in conjunction with products not recommended by Puregraft LLC.
  - Inspect the Puregraft Easel and Puregraft Slider before each use for sign of corrosion or deformities that would adversely affect its function or safety of use.
  - Inspect anti-slip pads on the base of the Puregraft Easel to ensure they are in place and intact.
- The physician is solely responsible, in all cases, when ordering use of this system. It is the physician's responsibility to use this system as intended.

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.

**European Reporting Requirements:** Any serious incident that has occurred in relation to the device should be reported to Bimini Health Tech and the competent authority of the Member State in which you, as the user and/or patient, are established.

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## EN: INSTRUCTIONS FOR USE

### DESCRIPTION

The Puregraft® Instrument Set (SST) was designed for use with the Puregraft 250 and 850 Systems, and is composed of the Puregraft Easel and Puregraft Slider. Both components are made of autoclavable stainless steel and can be sterilized as described in the procedures below. The Easel serves as a sterile stand to hold the Puregraft Bag upright and the Slider aids the user in extracting the graft from the Puregraft System once it's ready for re-injection. This is non-surgically invasive and intended to be cleaned and steam sterilized prior to use.

### CONTENTS

The Puregraft Instrument Set (SST) consists of the following:

DESCRIPTION	QUANTITY
Puregraft Easel (SST)	1 each
Puregraft Slider (SST)	1 each

### INSTRUCTIONS FOR USE

#### 1. CLEANING PROCEDURE (automated)

The recommended machine cleaning should follow a program similar to the Vario TD program from Miele, which has the following parameters:

- 2 minutes pre-cleaning with cold tap water
- Drain
- 5 minutes cleaning at 55°C with 0,5% solution of the
  - Drain
  - 3 minutes rinsing with cold demineralized water
  - Drain
  - 2 minutes final rinsing with cold demineralized water
  - Drain
  - A slightly alkaline cleaner like the neodisher MedClean (Dr. Weigert) is recommended.

#### 1.2 CLEANING PROCEDURE (manual)

- Prepare an enzymatic detergent according to the detergent manufacturer's recommendations.
- Fully immerse the device(s) in the prepared detergent and allow them to pre-soak for a minimum of 5 minutes.
- Thoroughly brush the device(s) ensuring all hard-to-reach areas are accessed. Use a syringe to aid in cleaning. Assume all visible soil is removed. If the detergent becomes grossly soiled, prepare a fresh batch as described in the first step.
- Thoroughly rinse the device(s) under running lukewarm tap water for a minimum of 1 minute. Use a syringe to aid in rinsing.

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