
# Puregraft Combined Adapter

## EN: INSTRUCTIONS FOR USE

### Indications For Use

The Reusable Puregraft ADAPT/PURE is intended to be used for the transfer of soft tissue and fluid between syringes and Puregraft products during general surgical and plastic surgery procedures. Reusable Puregraft ADAPT/PURE are intended for use in soft tissue surgical procedures and as advised by a qualified surgeon. This device is non-surgically invasive and intended to be cleaned and steam sterilized prior to use.

### Instructions For Use

Clean and sterilize all devices prior to use. Inspect the Puregraft ADAPT/PURE to ensure that there is no damage or wear. Connect the Puregraft ADAPT/PURE to the appropriate syringes as applicable.

### Contraindications

- Do not use Puregraft ADAPT/PURE if they do not meet the requirement(s) for use in the intended surgical procedure(s).
- Whenever and wherever hospital policy and procedure(s) take precedence over this protocol. After use in surgery, always clean the adapters prior to sterilization. Failure to use approved cleaning & sterilization procedures may void the Product Warranty.

### Warnings

- Reusable adapters are supplied non-sterile and require cleaning and sterilization prior to use. For cleaning and sterilization instructions, refer to this document.
- Do not use alkaline-based cleaning agents to clean Puregraft ADAPT/PURE. Alkaline-based (pH>7) detergents and solutions corrode metals. The center hole requires particular attention during cleaning.
- This document pertains to reusable devices only and these instructions are not intended for devices labeled "For Single-Use Only." Puregraft single-use devices must not be reused or reprocessed, as they are not designed to perform as intended after the first usage. Please refer to the device label and the package insert to identify the adapters as "For Single-Use Only" or "Reusable."
- Since reusable Puregraft ADAPT/PURE are supplied non-sterile and require cleaning and sterilization before use, the appropriate sections of these instructions may be applied unless other specific instructions are provided in the package insert. The healthcare facility should ensure that the combination of cleaning and sterilization used actually results in devices that are safe for use in surgery.
- Alternative methods of processing may be equally suitable. In the event that National cleaning and sterilization requirements conflict with these instructions, the National cleaning and sterilization requirement should be followed as applicable.

### Caution

Before using any Puregraft ADAPT/PURE, visually inspect it to ensure it is appropriate for the desired surgical application. Reusable connectors are vulnerable to repeated autoclaving (heat exposure to metals) that may cause it to wear out over time. Puregraft ADAPT/PURE removed from their mated syringes during surgery may

cause their mated syringes to produce syringe particulate. Before, during and after use, always

inspect the devices and syringes to ensure properly mated surfaces and fitness of use.

### Limitations and restrictions on reprocessing

Repeated heat exposure and reprocessing affects these adapters. Decommissioning of the devices is normally determined by wear and damage due to use.

### OPERATING INSTRUCTIONS

1. Connect a Toomey syringe or a catheter tip syringe to the larger end of the adapter.
2. Connect a Luer lock syringe to the smaller end of the adapter.
3. Transfer soft tissue and fluid between syringes as needed.
4. Puregraft Combined Syringe Adapter may also be used as a decant stand by placing the connected syringe upright.

### Instructions Point of Use

1. Remove excess blood and debris with disposable cloth and/or non-shedding wipe.
2. Puregraft ADAPT/PURE should be cleaned within 30 minutes of use to minimize the potential for drying prior to cleaning.
3. Disconnect adapters from their mated syringes prior to cleaning.
4. Prepare to pre-soak the adapters in a basin large enough to accommodate the adapters using an enzymatic cleaner according to the cleaning agent instructions.
- Preparation for Decontamination
5. Place devices in a solution as per enzyme solution manufacturer instructions using lukewarm (22-43°C) tap water.
6. Completely submerge all adapters in the enzyme solution and soak for at least 20 minutes.

### Cleaning: Manual

7. Remove the devices from the enzyme solution and rinse in lukewarm (22-43°C) tap water.
8. Clean the inside of the adapter surfaces with a soft-bristle brush and lukewarm (22-43°C), tap water mixed with the enzymatic cleaning agent as described by the manufacturer of the cleaning agent
9. Clean the inside lumens of the adapter using a soft-bristled, nylon brush until all visible debris has been removed. Using a syringe with a volume of 60ml, thoroughly and aggressively flush adapter lumens, holes and other difficult to reach areas.
10. Repeat at least 3 times with lukewarm (22-43°C) tap water-enzyme solution at a concentration recommended by the manufacturer.
11. Finish by flushing the adapter at least 3 times with a syringe (60ml) filled with ultrafiltered, RO, DI, distilled and/or demineralized water.
12. Place prepared enzymatic cleaning agent in an ultrasonic cleaner.
13. Completely submerge adapters in the cleaning solution and sonicate for 10 minutes at 40-45 kHz.
14. Rinse instrument in ultra-filtered, RO, DI, distilled and/or demineralized water for at least 3 minutes or until there is no sign of blood or soil on the device in the rinse stream. Repeat the sonication and rinse steps above.
15. Thoroughly dry moisture from the instrument with a clean, absorbent and non-shedding wipe.

### Cleaning: Automated

16. Remove the devices from the enzyme solution and rinse in lukewarm (22-43°C) tap water.
17. Transfer the adapter(s) into the automatic washer/disinfectant for processing.
18. Orient the adapter(s) to facilitate drainage in a way to allow flushing.
19. Select the cycle and ensure the following set of cycle parameters are properly programmed.
20. Remove the adapter(s) from the automatic washer/ disinfectant and dry using a clean soft non-shedding cloth and filtered pressurized air at < 40 psi.
21. Visually inspect the adapter(s) for cleanliness and repeat the cleaning procedure if necessary.

### Drying:

22. Wipe dry with a non-shedding wipe.
23. Ensure cleaned adapters are properly cleaned and dry prior to autoclaving by visually inspecting each instrument.

### Maintenance, Inspection and Testing

24. Inspect each instrument to ensure it is clean and has been sterilized before use.
25. Visually inspect each instrument for damage and signs of distortion.
26. If any sign of damage or distortion, call Puregraft immediately for repair or replacement.

### Packaging

Once cleaned and inspected, the dry Puregraft ADAPT/PURE should be double wrapped according to AAMI/ CSR technique. The packaging for terminally sterilized instruments should fulfill requirements EN ISO 11607.

### Sterilization

Steam sterilization (moist heat) is recommended.

### Pre-Vacuum

Wrapped Goods Cycle 4 minutes @ 132° C / 270° F. Dry Time: 20 to 30 minutes, according to load size.

### Storage

Identify and store sterilized adapters in accordance with EN ISO 11607.

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.

**European Reporting Requirements:** Any serious incident that has occurred in relation to the device should be reported to Bimini Health Tech and the competent authority of the Member State in which the you, as the user and/or patient, are established

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