Puregraft® 50

INSTRUCTIONS FOR USE
IN ALI COUNTRIES EXCEPT UNITED STATES AND CANADA
Puregraft® 50
EN: INSTRUCTIONS FOR USE

This product is certified as a medical device in the European Union under the Medical Device Regulation 2017/745 by SGS CE1639, exclusively for the indications of autologous fat transfer. Other non-medical uses as stated anywhere both containing ascribed to this device are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for these purposes.

INTENDED PURPOSE

The Puregraft® 50 is intended for use in the following surgical specialties when the aspiration of soft tissue is desired: plastic and reconstructive surgery, gastrointestinal and affiliated organ surgery, urological surgery, general surgery, orthopedic surgery, gynecological surgery, thoracic surgery, and laparoscopic surgery.

DESCRIPTION

The Puregraft® 50 is a sterile, single use, closed device intended for the preparation and delivery of autologous fat grafts back to the same patient for plastic and reconstructive surgery applications.

INTENDED USERS, PATIENTS, AND ENVIRONMENT

This device is intended to be used within a clinical setting by medical practitioners on patients who require the transfer of harvested autologous adipose tissue.

INDICATIONS FOR USE

The Puregraft 50 is indicated for use in the harvesting, filtering and transferring of autologous fat tissue for re-injecting back into the same patient for Autologous Fat Transfer (AFT) procedures. The Puregraft 50 is intended for use in the following surgical specialties when the aspiration of soft tissue is desired: plastic and reconstructive surgery, gastrointestinal and affiliated organ surgery, urological surgery, general surgery, orthopedic surgery, gynecological surgery, thoracic surgery, and laparoscopic surgery.

CONTRAINDICATIONS

The following are contraindications for use of the Puregraft 50:

- For autologous use only.
- Patients with the presence of any disease processes that adversely affect wound healing, and poor overall health status of the individual.

STERILITY

The Puregraft 50 is sterilized with gamma irradiation.

STORAGE AND HANDLING

Use only if the sterilization indicator is red.

WARNINGS

- Single Use Only – Do Not Re-use.
- For Autologous Use Only.
- The surgeon should be familiar with the operation of the Puregraft 50 and the surgical procedure prior to performing the surgery.
- Do not re-use the Puregraft 50 during the same procedure. Re-use may compromise the sterility and adversely affect performance of the device.
- Re-use of a single-use device presents significant risk to the proper operation of the device due to absence of the manufacturer's Instructions For Use (IFU).
- Re-use or cleaning of a single-use device may result in harm to the patient due to unknown chemical interactions with the device and/or exposure of residual cleaning agents to the patient.
- Re-use of a single-use device poses serious risks for the patient to include loss of sterility which could lead to infection.
- Re-use of a single-use device poses serious risks for the patient to include the transmission of infectious disease from one patient to another.
- In addition to infectious disease transmission risks, the re-use of a single-use devices poses serious risks for the patient to include immediate and/or premature mechanical failure of the device due to material degradation which could lead to reduced or impaired performance of the product up to and including catastrophic failure of the device.
- Re-sterilization of this device may result in significant damage or loss of integrity to the device.
- Do not overfill the Puregraft 50 with more than 100 mL of total volume of tissue (50 mL) and/or washing solution (50 mL).
- Immediately stop the introduction of any fluid or tissue into the Puregraft 50 if the material of the bag begins to form leading to loss of integrity of the bag and the creation of a biohazard to the operating staff.

PRECAUTIONS

- Inspect the Puregraft 50 package for signs of damage prior to using.
- The Puregraft 50 is single for use only and must be used immediately after removal from the sterile packaging. Discard the Puregraft 50 if the packaging is damaged, opened or the sterilization indicator is invalid.
- Aseptic techniques must be used to minimize the possibility of contamination of the Puregraft 50.
- Treat all blood and fluids using Universal Precautions (Blood-borne Pathogen Precautions).
- Over inflating or over filling of the Puregraft 50 may result in a biohazard to the operating staff and/or device malfunction.
- Adipose tissue harvested using cannulae larger than 4 mm may cause clogging.
- Extracting graft too quickly may cause the Puregraft 50 to collapse.
- Potential adverse events found in the literature for fat transfer include abscess, bleeding/ Blood loss, Cellulitis, Cysts, Fat necrosis, graft loss, granuloma, hematoma, infection, nodules/palpable mass, seroma
- Potential device malfunctions found in the literature for fat transfer include blockage; burst container or vessel; collapse; contamination; detachment of device or device component; device markings/labelling problem; blood/blood leak, hole, tear, or rip in the material; hole, tear, or rip in the barrier esthetic; leak; loose or intermittent connection; material separation; malfunction causing all the fat to be evacuated into the non-sterile waste collection canister; device malfunctions; presence of bacterial endotoxins; pressure problem; and splash.
- The following hazardous situations were identified as having a high severity via risk analysis:
 - Sterilization - use of a product which is not sterile could lead to infection or death. The Puregraft 50 device is sterilized to a Sterility Assurance Level of 10⁻⁶.
 - Microbial Contamination - use of a product which has been contaminated and as compromised sterility could lead to infection. The Puregraft 50 device could be packaged using industry-standard Typey packaging, and testing has been performed to verify that the product packaging does not degrade after worst-case sterilization conditions and aging.
 - Biocompatibility - use of a materials which are not biocompatible could result in an adverse reaction, leading to irritation, injury, sickness, or disease. Bimini Health Tech performs Biocompatibility testing to demonstrate biocompatibility of the materials used in the Puregraft 50 device.
 - Graft rejection - use of allogenic or xenogenic tissue could result in graft rejection. The Puregraft 50 device is labeled for autologous fat grafting only.
 - Misuse - misuse of the Puregraft 50 device could include intravenous applications, allogenic processing, use of non-adipose tissue, and xenotransplantation. Misuse is warned against in this Instructions for Use.

OPERATING INSTRUCTIONS

MATERIALS REQUIRED FOR PROCESSING	QUANTITY
Puregraft 50	1
Minimum 100 mL Lactated Ringer's (not provided)	1
100 mL syring(s) (not provided)	
10mL or smaller	2 or more
60 mL or smaller	2 or more

Note: Follow standard surgical sterility procedures when processing with the Puregraft 50.

Note: Refer to step symbols **[1]****[2]****[3]****[4]** when following operating instructions.

- PREPARATION:
- Open the Tyvek® pouch containing the sterile Puregraft 50.

TISSUE ADJUSTMENT

- Introduce 30-50 mL of harvested tissue into the Puregraft 50 through the "Tissue/Wash" port using a luer lock syringe.

GRAFT WASHING AND PROCESSING

- FIRST WASH:
- Aspirate 50 mL of Lactated Ringer's to the Puregraft 50 through the "Tissue/Wash" port, regardless of input tissue volume.
 - Agitate contents in the Puregraft 50 by rocking the device back and forth for a minimum of 15 seconds.
 - Place the Puregraft 50 on a flat surface from side up for a minimum of 1 minute to allow contents to separate.
 - Attach an empty luer lock syringe to the "Drain" port.
 - Hold bag at an upright, angled position and slowly extract all waste fluid through the "Drain" port.

SECOND WASH:

- Repeat steps 3 – 7. Two washes are required for optimal graft output.

GRAFT REMOVAL

- Extract and discard any excess fluid trapped in the "Tissue/Wash" port.
- Hold slider at the sides and pull down to direct the graft towards the ports for extraction.
- Extract the graft through the "Tissue/Wash" port using a new luer lock syringe.
- Once the graft is removed from the Puregraft 50, properly dispose all components using Universal Precautions for Blood-borne Pathogens.

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PERFORMANCE CHARACTERISTICS

- 10-25% of tissue output given 30-35 mL of harvested tissue.
- Withstands up to 3.0 PSI for up to 10 seconds.
- Operation between 15 °C – 30 °C.

European Reporting Requirements: Any serious incident that has occurred in relation to the use of the device should be reported to Bimini Health Tech and the competent authority of the Member State in which the you, as the user and/or patient, are established.

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